



## Certificate of Analysis

Page 1 of 1

<b>Client:</b>	Pomona Group Limited	<b>Lab No:</b>	3177063	SPv1
<b>Contact:</b>	Ian Snelgrove	<b>Date Received:</b>	18-Feb-2023	
	C/- Pomona Group Limited	<b>Date Reported:</b>	23-Feb-2023	
	PO Box 21806	<b>Quote No:</b>	122346	
	Henderson	<b>Order No:</b>		
	Auckland 0650	<b>Client Reference:</b>		
		<b>Submitted By:</b>	Ian Snelgrove	

### Sample Type: Fomite

<b>Sample Name:</b>	Latex 330S	Latex 330M	Latex 330L	Latex 330XL
<b>Lab Number:</b>	3177063.1	3177063.2	3177063.3	3177063.4
Total Coliforms	cfu / item	< 1 #1	< 1 #1	< 1 #1

### Analyst's Comments

#1 Statistically estimated count based on the theoretical countable range for the stated method.

## Summary of Methods

The following table(s) gives a brief description of the methods used to conduct the analyses for this job. The detection limits given below are those attainable in a relatively simple matrix. Detection limits may be higher for individual samples should insufficient sample be available, or if the matrix requires that dilutions be performed during analysis. A detection limit range indicates the lowest and highest detection limits in the associated suite of analytes. A full listing of compounds and detection limits are available from the laboratory upon request. Unless otherwise indicated, analyses were performed at Hill Laboratories, 28 Duke Street, Frankton, Hamilton 3204.

### Sample Type: Fomite

Test	Method Description	Default Detection Limit	Sample No
Total Coliforms	Sample item/s rinsed in sterile buffer, Membrane Filtration, Resus on TSA, Count on mEndo Les agar, Incubation at 35°C for 20.5-22 hours, Confirmation. APHA 3.7 5 <sup>th</sup> Ed (modified).	1 cfu / item	1-4

These samples were collected by yourselves (or your agent) and analysed as received at the laboratory.

Testing was completed on 23-Feb-2023. For completion dates of individual analyses please contact the laboratory.

Samples are held at the laboratory after reporting for a length of time based on the stability of the samples and analytes being tested (considering any preservation used), and the storage space available. Once the storage period is completed, the samples are discarded unless otherwise agreed with the customer. Extended storage times may incur additional charges.

This certificate of analysis must not be reproduced, except in full, without the written consent of the signatory.

Shobhna Ram BSc  
Operations Support - Microbiology